



L100R

Keep refrigerated

Perforated card design is the subject of a pending patent application

Safety information

Slightly hazardous (irritant, sensitizer) in case of skin and/or eye contact, always wear gloves and safety glasses.

Description:

LyseNow® perforated Cards are thick filter paper card that were treated with proprietary chemical formula. The card preserves DNA/RNA integrity at ambient temperature. Each card has seven perforated 3 mm discs for easy detach with pipette tips. Each card is labeled with individual bar code for easy logging and tracking.

Kit contents

item	description	quantity
LyseNow® Perforated Card	Individually packaged in zip bag	25

Protocol

- Sample application on card**
 - Directly drop up to 100uL of fluid on the center of perforated area;
 - Or, collect sample on swab, press and roll swab onto the perforated discs;
- Dry the card** on a portable Card Drying station (Cat. # U100) for about 30 minutes or at ambient temperature for about two hours.
- Card storage**
 - Return the card to its original zip bag.
 - DNA are stable at ambient temperature for at least a year;
 - RNA are stable at ambient temperature for at least a week;
 - If accessible, store cards at 4 °C or -20 °C for longer storage.
- Disc translocation**
 - Push out seven discs into an eppendorf tube using a sterile pipette tip;

5. RNA recovery


- Clean Sample RNA (oral/nasal fluid, serum, etc)**
 - Submerge seven discs in 100uL water*, vortex 3 times, each for 5 sec;
 - Use 100uL of 5% Chelex 100 for better RNA recovery
 - Spin the tube at top speed for 1 min;
 - Transfer the supernatant containing recovered RNA to a new eppendorf tube.
- Tough Sample RNA (blood, feces, soil, etc)**
 - Submerge seven discs in 350uL Trizol, or phenol:chloroform, or lysis buffer from RNA purification kits;
 - Vortex the tube at top speed for 5 minutes;
 - Transfer the supernatant to a new eppendorf tube;
 - Follow the protocol of selected RNA purification method to further purify RNA.

❖ See back for continue

6. DNA recovery

- Clean Sample DNA (oral/nasal fluid, serum, etc)**
 - Submerge seven discs in 100uL water, heat to 95C for 30 min;
 - Spin the tube at top speed for 1 min;
 - Transfer the supernatant containing recovered DNA to a new eppendorf tube.
- Tough Sample DNA (blood, feces, soil, etc)**
 - Submerge seven discs in 500uL water, vortex at top speed for 3x5 sec, discard supernatant;
 - Add 100uL of nuclease free water;
 - Heat the tube in a 95 °C heating block for 30 min;
 - Vortex the tube at top speed for 5 sec;
 - Spin the tube at top speed for 1 min;
 - Transfer the supernatant containing recovered DNA to a new eppendorf tube.

Brief workflow



1

2

3

4

5