



S200R

Keep refrigerated

Safety information

Slightly hazardous (irritant, sensitizer) in case of skin and/or eye contact, always wear gloves and safety glasses.

Description:

RNA*Sound*™ RNA Sampling Card unprecedentedly stabilizes RNA on filter paper with the proprietary impregnating lysis buffer. It features:

- Room temperature RNA sample collection, storage and transportation;
- Immediate inactivation and thus safe handling of infectious agents.

ReadyPunched™ format (Patent pending) eliminates the tedious card punching and the risk of cross contamination.

Kit contents

item	description	quantity
RNA <i>Sound</i> ™ ReadyPunched™ RNA Sampling Cards	Individually packaged in zip bag with desiccant	25

Protocol

- Sample preparation**
 - 1) Fluid sample (blood, serum, saliva, culture supernatant) can be applied directly on card;
 - 2) Feces or dirt: Add 0.5 g sample to 1 mL PBS and vortex for 3 min; centrifuge at 100 x g for 1 min and apply supernatant to the card.
 - 3) Cells or bacteria cultures:
 - (For adherent cells) Detach cells and inactivate trypsin;
 - Cells pelleted down;
 - Cells washed with 1XPBS;
 - Cells resuspended in 1XPBS

1

2

3

4


2. Sample application on card

- 1) Directly drop sample on the perforated discs on the card;
 - The two perforated discs on each card are for a same sample.
- 2) Or, collect sample on a swab, and press and roll the contents of the swab onto the perforated discs on the card;
 - The stability of sample RNA is not guaranteed outside the two perforated discs.
- 3) Dry the card on a portable Card Drying station (Cat. # U100) for about 10 minutes or at room temperature for about an hour.
- 4) Return the card to its original zip bag with desiccant.
 - RNA are stable at room temp for at least one week;
 - If accessible, store cards at 4 °C or lower for longer storage.

3. RNA elution

- 1) Take the card out of the zip bag;
- 2) Further dry the card if necessary on portable Card Drying station (Cat. # U100), or in air for an hour;
- 3) Push out the two perforated discs into a 1.5 mL eppendorf tube using a sterile pipette tip;
 - If the disc hangs on the card, push the disc against the tube wall, and pull the card to detach the disc.
- 4) submerge the paper pieces with standard volume of Trizol, or phenol:chloroform, or lysis buffer from RNA purification kits;
- 5) Hand hold the tube and vortex at top speed for 5 minutes, or shake multiple tubes in an Eppendorf Thermomixer at top speed and 60C for 2 min.
- 6) Transfer the elute to a new eppendorf tube;
- 7) Follow the protocol of selected RNA purification method to purify RNA.
- 6) Use elute for immediate RT-PCR or store at -20C for future use.

Brief workflow



- Squeeze and roll the cotton tip on perforated discs;
- Dry the card on portable Card Drying station (Cat. # U100) for 10 minutes; or in air for 1 hour

- RNA are stable for one week;
- Push out discs to a 1.5 mL eppendorf tube

Add ~350 μL of RNA kit lysis buffer, vortex for 5 minutes to elute RNA, follow the protocol of selected RNA purification method to purify RNA.

5

6

7